

REFRACTOMETRY ISO SOLUTION °BX 10,0 / nD 1.34782 PACK 15x4 ml vials

Art. ID	CF-CHFISOBX1015F
Unit	15 x 4 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Refractometry standard solution (°BX / nD) with ISO traceable certificate /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 20 °C /// We can produce CUSTOMIZED SOLUTIONS at any point between 0 and 60 Bx (between 1.33299 and 1.44193nD) (quotation on request).