

**REFRACTOMETRY ISO SOLUTION °BX 12,00 / nD 1.35093 PACK 15x4 ml vials**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | CF-CHFISOBX1215F                  |
| Unit            | 15 x 4 mL                         |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Refractometry standard solution (°BX / nD) with ISO traceable certificate /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 20 °C /// We can produce CUSTOMIZED SOLUTIONS at any point between 0 and 60 Bx (between 1.33299 and 1.44193nD) (quotation on request).