

**REFRACTOMETRY ISO SOLUTION °BX 20,00 / nD 1.36384 PACK 15x4 ml vials**

Art. ID	CF-CHFISOBX2015F
Unit	15 x 4 mL
Deliverydetails	No Dangerous Good /not restricted

**Description**

Refractometry standard solution (°BX / nD) with ISO traceable certificate /// Product shelf life minimum 1 year from purchase /// The values given refer to the reference temperature of 20 °C /// We can produce CUSTOMIZED SOLUTIONS at any point between 0 and 60 Bx (between 1.33299 and 1.44193nD) (quotation on request).