

**REFRACTOMETRY ISO SOLUTION °BX 30,00 / nD 1.38115 250 ml bottle**

Art. ID	CF-CHFISOBX30250
Unit	250 mL
Deliverydetails	No Dangerous Good /not restricted

**Description**

Refractometry standard solution (°BX / nD) with ISO traceable certificate /// Product shelf life 6 months from purchase /// The values given refer to the reference temperature of 20 °C /// We can produce CUSTOMIZED SOLUTIONS at any point between 0 and 60 Bx (between 1.33299 and 1.44193nD) (quotation on request).