

### **Custom Solution - 4 components in HNO3 5%**

Art. ID LM24-CUS-71778  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/l		
	Boron (B)	[7440-42-8]	100	mg/l		
	Barium (Ba)	[7440-39-3]	100	mg/l		
	Silicon (Si)	[7440-21-3]	500	mg/l		