

Aripiprazole Acid Analog (4-((2-oxo-1,2,3,4-tetrahydroquinolin-7-yl)oxy)butanoic acid)

Art. ID USP-1A02470
Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-((2-oxo-1,2,3,4-tetrahydroquinolin-7-yl)oxy)butanoic acid	[58899-27-7]				