

**\*\*\*closed\*\*\*Round robin on recycling material (proficiency test)**

Art. ID	L-RR-2204RC
Unit	Two different pieces of Aluminium alloy profile-scrap
Deliverydetails	No Dangerous Good /not restricted

**Description**

The round robin test will be carried out on two different pieces of Aluminium alloy profile-scrap. Both alloys are within the same alloy group. Purpose of this round robin is the identification of the alloy type as well as the chemical composition. Furthermore, it shall give the participants information about if their analysis can differentiate between the two alloys ///The material shall be prepared and analyzed by the common method of the reception control. For example LIBS, XRF, OES or ICP. The alloy number and the method used shall be reported additionally to the results ///After the completion of the test, each participant will receive a copy of the test's evaluation and a certificate of participation, according to the DIN EN 38402-42 standard. An evaluation on the identification and the used methods will be done as well /// Due to the nature of the material and the use of different analytical methods, the round robin test is only carried out after a certain minimum number of participants has been reached in order to achieve statistical robustness.