

Steel - Stainless, CRM, 309, Chips

Art. ID IARM-Fe309-18-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,0065	%		
	Carbon (C)	[7440-44-0]	0,066	%		
	Cobalt (Co)	[7440-48-4]	0,248	%		
	Chromium (Cr)	[7440-47-3]	22,42	%		
	Copper (Cu)	[7440-50-8]	0,43	%		
	Manganese (Mn)	[7439-96-5]	1,61	%		
	Molybdenum (Mo)	[7439-98-7]	0,357	%		
	Niobium (Nb)	[7440-03-1]	0,021	%		
	Nickel (Ni)	[7440-02-0]	12,2	%		
	Phosphorus (P)	[7723-14-0]	0,029	%		
	Sulfur (S)	[7704-34-9]	~0,002	%		
	Antimony (Sb)	[7440-36-0]	0,0018	%		
	Silicon (Si)	[7440-21-3]	0,3	%		
	Tin (Sn)	[7440-31-5]	0,012	%		
	Vanadium (V)	[7440-62-2]	0,073	%		
	Tungsten (W)	[7440-33-7]	0,063	%		