

Aripiprazole 4,4'-dimer
(7,7'-((((ethane-1,1-diylbis(2,3-dichloro-4,1-phenylene)))bis(piperazine-4,1-diyl)))bis(butane-4,1-diyl)))bis(oxy))bis(3,4-dihydroquinolin-2(1H)-one)

Art. ID USP-1A02460

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aripiprazole 4,4'-dimer	[1797986-18-5]				