

Cast Iron - High Alloy, CRM, last of stock, Block ~36 mm x ~36 mm x ~18 mm

Art. ID VS-ChG 44
Unit Block
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Carbon (C) | [7440-44-0] | 1,24 | % | | |
| | Chromium (Cr) | [7440-47-3] | 25,44 | % | | |
| | Copper (Cu) | [7440-50-8] | 2,27 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,87 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,035 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,175 | % | | |
| | Phosphorus (P) | [7723-14-0] | ~1,2 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,076 | % | | |
| | Silicon (Si) | [7440-21-3] | 1,5 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,104 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,079 | % | | |