

Custom Solution - 5 components 100 mg/L each of in HNO3 5%

Art. ID LM24-CUS-71790
Unit 2 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Zirconium (Zr)	[7440-67-7]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		