

## **Custom Solution - KIT of 27 components 100 µg/mL each of in Methanol**

Art. ID LM24-CUS-71795  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fenoprop	[93-72-1]	100	µg/mL		
	2,4,5-T	[93-76-5]	100	µg/mL		
	Aldrin	[309-00-2]	100	µg/mL		
	Chlordane	[57-74-9]	100	µg/mL		
	4,4'-DDD	[72-54-8]	100	µg/mL		
	4,4'-DDE	[72-55-9]	100	µg/mL		
	4,4'-DDT	[50-29-3]	100	µg/mL		
	Dieldrin	[60-57-1]	100	µg/mL		
	Endrin	[72-20-8]	100	µg/mL		
	Heptachlor-exo-epoxide (cis-, isomer B)	[1024-57-3]	100	µg/mL		
	Hexachlorobenzene	[118-74-1]	100	µg/mL		
	Pentachlorobenzene	[608-93-5]	100	µg/mL		
	1,2,4,5-Tetrachlorobenzene	[95-94-3]	100	µg/mL		
	alpha-HCH	[319-84-6]	100	µg/mL		
	beta-HCH	[319-85-7]	100	µg/mL		
	delta-HCH	[319-86-8]	100	µg/mL		
	Isodrin	[465-73-6]	100	µg/mL		
	Kelevan	[4234-79-1]	100	µg/mL		
	Chlordecone-hydrate	[143-50-0]	100	µg/mL		
	gamma-HCH (Lindane)	[58-89-9]	100	µg/mL		
	Methoxychlor	[72-43-5]	100	µg/mL		
	Mirex	[2385-85-5]	100	µg/mL		
	Perthane	[72-56-0]	100	µg/mL		
	Quintozene	[82-68-8]	100	µg/mL		
	Strobane	[8001-50-1]	100	µg/mL		
	Isobenzan (Telodrin)	[297-78-9]	100	µg/mL		
	Campechlor (Toxaphene)	[8001-35-2]	100	µg/mL		

