

**Customized Set of SURAGUS certified calibration reference samples for 0,03 - 1000 Ohm/sq (coated PET) (SURAGUS manufacturer certification) (on FR-4 Substrate)**

Art. ID	SUR-1000XXX-SET-FR-4
Unit	Customized (Size: from 1 inch (round) and from 10 x 10 mm (rectangular))
Deliverydetails	No Dangerous Good /not restricted

**Description**

Calibration of electrical characterization systems / Validation of system accuracy to meet ISO 9001 and ISO/IEC 17026