

Assy for incoming gases for Autofluxer of mowitech GmbH includes following parts: pre mounted on a new bracket, pre installed and tested, includes also the new cable for the three valves

Art. ID MOWI-SET-ASSY-FOR-INCOMING-GASES
Unit ea
Delivery details No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
MOWI-LTM AFX-VO2-ASS						
Valve O2- Assy with met						
al tube for Autofluxer						
of mowitech GmbH						
MOWI-LTM AFX-VG2-ASS						
SW: Valve G2- Assy leak						
tested, with connector						
s and gas check valve f						
or Autofluxer of mowite						
ch GmbH						
MOWI-LTM AFX-VL10-ASS						
Valve pneumatic cylind						
er with connectors for						
Autofluxer of mowitech						
GmbH						
MOWI-LTM-tube-con-(Air)						
: Tube connector air, o						
Id style for Autofluxer						
of mowitech GmbH						
MOWI-LTM-tube-con-(Gas)						
: Tube connector for ga						
s, old style for Autofl						
uxer of mowitech GmbH						
MOWI-LTM-tube-con-(O2):						
d=6mm - Tube connector						
for oxygen for Autoflu						
xer of mowitech GmbH						
MOWI-LTM-tube-con-(quic						
k-air): Tube connector						
air quick connect for A						

utofluxer of mowitech G

mbH

MOWI-LTM AFX-V-1/4-4: T

ube union d=4mm, G1/4 s

traight (air) for Autof

luxer of mowitech GmbH

MOWI-LTM-tube-con-(gas-

1/2"): Tube connector

for gas 1/2" to 8mm fo

r Autofluxer of mowitec

h GmbH

MOWI-LTM AFX-cylinder:

Pneumatic cylinder asse

mbly for Autofluxer of

mowitech GmbH

MOWI-LTM AFX-Air-regula

tor: Internal air regul

ator with tube connecto

rs and gauge for Autofl

uxer of mowitech GmbH