

Custom Solution - Heptachlor 200 µg/mL in Iso-octane

Art. ID LM24-CUS-71005-1x10ML
Unit 1 x 10 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Heptachlor	[76-44-8]	200	µg/mL		