

Custom Solution - 11 components in HNO3 5%

Art. ID LM24-CUS-71848
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Nickel (Ni)	[7440-02-0]	50	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Barium (Ba)	[7440-39-3]	200	mg/L		
	Strontium (Sr)	[7440-24-6]	200	mg/L		
	Zinc (Zn)	[7440-66-6]	200	mg/L		
	Phosphorus (P)	[7723-14-0]	200	mg/L		
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Lithium (Li)	[7439-93-2]	1000	mg/L		