

## **Custom Solution - 14 components in HNO3 5% + Extended Shelf-Life (Option D): + 8 months**

Art. ID LM24-CUS-71857  
Unit 1 x 250 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Conservation 8 months + Extended Shelf-Life (Option D) + 8 months. (4 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	40	mg/L		
	Arsenic (As)	[7440-38-2]	4	mg/L		
	Cadmium (Cd)	[7440-43-9]	0,4	mg/L		
	Chromium (Cr)	[7440-47-3]	4	mg/L		
	Copper (Cu)	[7440-50-8]	400	mg/L		
	Iron (Fe)	[7439-89-6]	80	mg/L		
	Manganese (Mn)	[7439-96-5]	40	mg/L		
	Nickel (Ni)	[7440-02-0]	8	mg/L		
	Phosphorus (P)	[7723-14-0]	8	mg/L		
	Lead (Pb)	[7439-92-1]	4	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Selenium (Se)	[7782-49-2]	10	mg/L		
	Silicon (Si)	[7440-21-3]	200	mg/L		
	Zinc (Zn)	[7440-66-6]	80	mg/L		