

**Custom Solution - 3 components 5 mg/L each of in HNO3 2% + Extended Shelf-Life (Option D): + 12 months**

Art. ID LM24-CUS-71861  
Unit 2 x 250 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months + Extended Shelf-Life (Option D) + 12 months. (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	5	mg/L		
	Arsenic (As)	[7440-38-2]	5	mg/L		
	Iron (Fe)	[7439-89-6]	5	mg/L		