

Custom Solution - 5 components 100 µg/mL each of in 2,2,4-Trimethylpentane (Isooctane)

Art. ID LM24-CUS-71864
Unit 3 x 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Decabromodiphenyl ether (PBDE 209)	[1163-19-5]	100	µg/mL		
	Tris(4-isopropylphenyl) Phosphate	[2502-15-0]	100	µg/mL		
	2,4,6-Tri-tert-butylphe nol	[732-26-3]	100	µg/mL		
	Hexachloro-1,3-butadien e	[87-68-3]	100	µg/mL		
	Pentachlorothiophenol	[133-49-3]	100	µg/mL		