

Custom Solution - Hg 3 mg/L in HNO3 5%

Art. ID LM24-CUS-71870
Unit 1 x 50 mL bottle
Deliverydetails No Dangerous Good UN3264 II E2

Description

Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Mercury (Hg) | [7439-97-6] | 3 | mg/L | | |