

Custom Solution - 16 components in HNO3 2%

Art. ID LM24-CUS-71891
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months. (6 months after opening).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 10 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 1 | mg/L | | |
| | Boron (B) | [7440-42-8] | 10 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 1 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 1 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 10 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 10 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 10 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 10 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 1 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 1 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 1 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 1 | mg/L | | |
| | Thallium (Tl) | [7440-28-0] | 1 | mg/L | | |
| | Uranium (U) | [7440-61-1] | 1 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 10 | mg/L | | |