

Custom Solution - 5 components 1000 mg/L each of in HNO3 5%, HF traces

Art. ID LM24-CUS-71221-2-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	1000	mg/L		
	Beryllium (Be)	[7440-41-7]	1000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Titanium (Ti)	[7440-32-6]	1000	mg/L		