

ReadyTube(R) Yeast extract agar For the enumeration of culturable microorganism in water samples in accordance with ISO 6222, pkg of (100 mL in 125 mL bottle with blue screw cap and 3 loci (10 bottles per box))

Art. ID SAF-1463520010

Unit EA

Description

ReadyTube(R) is ready to use media to provide the highest level of quality and testing confidence for improving accuracy and significantly reducing the risk of false results. ReadyTube(R) Yeast extract agar is a highly nutritive medium recommended for the plate count of microorganisms in water. The media contains Yeast extract and peptone to provide nitrogenous compounds, vitamin B complex, and other growth nutrients. The Yeast extract agar formulation is in accordance to ISO 6222, the test strains according to ISO 11133 and released through an ISO 17025 accredited laboratory.