

OmniPur® TES TES, CAS 7365-44-8 , is a Zwitterionic buffer that is useful in the pH range of 6.8-8.2.

Art. ID SAF-5051470001
Unit EA
Deliverydetails No Dangerous Good /not restricted

Description

A zwitterionic buffer useful in the pH range of 6.8-8.2. Has a pKa 7.40 at 25C., A zwitterionic buffer (Absorbance: <=0.12 at 260 nm, 1 M, H2O). Has a pKa of 7.5 at 20C. Useful in the pH range of 6.8-8.2.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	TES	[7365-44-8]				