

2,3-Epoxypropyl methacrylate Technipur(TM), for synthesis

Art. ID SAF-8006092500
Unit EA
Deliverydetails UN2922 (6.1) III E1

Description

2,3-Epoxypropyl methacrylate is a methacrylic acid/glycidol ester. It is used as a monomer to produce polymers based on glycidyl methacrylate (GMA). It is also a dual functionality (epoxy and methacrylate) monomer, thus offering the opportunity for chemical modification of the parent polymer for various applications.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Glycidyl methacrylate	[106-91-2]				