

Rink Amide AM resin (100-200 mesh) Novabiochem(R)

Art. ID SAF-8551300100

Unit EA

Description

This support comprises the modified Rink amide linker attached to aminomethylpolystyrene, and is an ideal tool for the Fmoc SPPS of peptide amides. Cleavage from this resin can be effected by a single step treatment with 95% TFA, providing peptide amides in high yields and purities.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-(2',4'-Dimethoxypheny	[183599-10-2]				
	I-Fmoc-aminomethyl)-phe					
	noxyacetamido-aminometh					
	yl resin					