

ESGRO Complete PLUS Clonal Grade Medium, 100ml The ESGRO Complete PLUS Clonal Grade Medium is a defined serum-free medium provided with a selective GSK3beta inhibitor to enhance viability of mouse embryonic stem cells & increased maintenance of pluripotency in the absence of serum & feeder cells.

Art. ID SAF-SF001-500P
Unit 1 x 500P

Description

ESGRO Complete PLUS Clonal Grade Medium is a defined serum-free medium provided with a selective GSK3beta inhibitor to enhance viability of mouse embryonic stem (ES) cells and increased maintenance of pluripotency in the absence of serum and feeder cells. This medium is intended for mouse ES cell expansion and growth at clonal density in serum-free conditions. The mouse ES cells should maintain germ line competency.