

MILLIPLEX(R) Mouse TH17 Magnetic Bead Panel - Immunology Multiplex Assay
Simultaneously analyze multiple Th17 cytokine and chemokine biomarkers with the Th17
Bead-Based Multiplex Assays using the Luminex technology, in human serum, plasma
and cell culture samples.

Art. ID SAF-MTH17MAG-47K-14C
Unit EA

Description

The MILLIPLEX(R) Mouse Th17 Bead panel enables you to focus on the therapeutic potential of the modulation of CD4 T-helper cell cytokines. Coupled with the Luminex(R) xMAP(R) platform in a bead format, you receive the advantage of ideal speed and sensitivity, allowing quantitative multiplex detection of dozens of analytes simultaneously, which can dramatically improve productivity. Panel Type: Cytokines/Chemokines