

Immobilon(R)-E PVDF Membrane 1 roll, 26.5 cm x 1.875 m, The first PVDF transfer Membrane for western blotting that wets out in transfer buffer or water. For chemiluminescence or chromogenic detection.

Art. ID SAF-IEVH00005

Unit EA

Description

Primary Binding Mechanisms: Electrostatic, hydrophobic, Filter Code: IEVH, Immobilon-E is the first PVDF membrane for Western blotting that can be wetted with aqueous transfer buffer. This eliminates the alcohol pre-wet step required with other PVDF membranes.