

Millipore Glass Microanalysis Filter Holder - Kit 25 mm, Glass frit membrane support, 15mL Funnel

Art. ID SAF-XX1012500
Unit EA

Description

Millipore Glass Microanalysis Filter Holders - Kit 25 mm, Glass frit membrane support, 15mL Funnel
The Millipore Glass Microanalysis filter holder is designed for vacuum filtration of small-volume liquid samples for particulate or biological contamination analysis. This item is a Kit composed of the following items:
(Catalog Number - Description)
XX1012514 - Millipore Glass funnel, 15 mL, 25 mm, Borosilicate, Ground glass seal
XX1012503 - Millipore Spring Clamp, 25 mm, Aluminum
XX1012502 - Millipore Glass Base and Stopper, 25 mm, Glass frit support
ALIGNED25 - Millipore Alignment Guide, 25 mm, 5/PK
Features & Benefits: . Volume graduations clearly indicate sample volume. Filter support available in fritted glass or stainless steel
Applications: Contamination Analysis (not for use with thick prefilters or flammable liquids)