

ech2o(R) ASM Lamp Protects stored water against biofilm buildup in Milli-Q(R) IQ/IX storage tanks

Art. ID	SAF-ASMUVLPA1
Unit	EA
Deliverydetails	No Dangerous Good

Description

The ech2o(R) ASM lamp protects the integrity of the stored water with regular exposure to bactericidal UV light. This ech2o(R) ASM lamp is a mercury-free UV LED (265 nm) that prevents the bacterial growth and biofilm formation on the inner surface of Milli-Q(R) IQ/IX storage tanks.