

**SMC(R) Human IL-1alpha High Sensitivity Kit kit sufficient for 96 wells, input: plasma, serum**

Art. ID SAF-03-0167-00

Unit EA

**Description**

The SMC(R) human IL-1alpha High Sensitivity Immunoassay Kit contains all reagents required to perform the assay. This assay takes advantage of a standard immunoassay workflow configured in a standard plate format. Assay protocols are similar to existing sandwich ELISA methods with two key differences: 1) Elution buffer disrupts the sandwich, separating the labeled detection antibody for quantification. 2) The Erenna and SMCxPRO(R) immunoassay systems detect analytes using Single Molecule Counting (SMC(R)) technology.