

MilliSentials Aliquoting Pipette Controller

Art. ID SAF-MISEPTC0N

Unit EA

Description

The MilliSentials Aliquoting Pipette Controller combines the automatic aliquoting function of a repeat pipettor with the basic functions of a standard pipette controller in a single device. The aliquoting function allows users to set an aliquot volume between 0.5-5.0 mL and automatically dispense accurate aliquots with the touch of a button. The device also provides touch sensitive aspirate and dispense buttons for standard pipette controller function and is compatible with all major brands of serological pipettes with volumes between 1-100mL. The MilliSentials Aliquoting Pipette Controller combines repeat pipetting and standard pipetting into one easy to use instrument. Accessories Included-. Wall Mount. 0.22 µm replaceable nozzle filter. Power cord