

Milliflex(R) Oasis Funnel Individually Packed, for use with Milliflex Oasis(R) systems, gridded white mixed cellulose esters (MCE) filter membrane, pore size 0.45 µm

Art. ID SAF-MMP000489

Unit EA

Description

Milliflex(R) Oasis Funnel is a pre-sterilized, ready-to-use filtration unit. The unit is comprised of: 100ml funnel, White gridded MCE filter membrane. The sealed integral membrane filter eliminates the need for membrane manipulation and handling, thus minimizing a major source of contamination and thereby improves bioburden testing reliability. The sealed membrane prevents the bypass phenomenon often found with other filtration systems. The Milliflex(R) filtration possesses almost twice the surface area of standard membranes enabling faster filtration and the funnel can be easily detached from the membrane without the use of any tools. The convex shape of the membrane after filtration demonstrates the integrity of the membrane and prevents air bubbles from being trapped when placed on culture media.