

**Phosphatase Inhibitor Cocktail Set V, 50X, Lyophilized The Phosphatase Inhibitor Cocktail Set V, 50X, Lyophilized controls the activity of Phosphatase. This small molecule/inhibitor is primarily used for Phosphorylation & Dephosphorylation applications.**

Art. ID	SAF-524632-1SET
Unit	1 x 1 SET
Deliverydetails	No Dangerous Good

**Description**

A cocktail of four phosphatase inhibitors for broad-spectrum inhibition of acid, alkaline, serine/threonine and protein tyrosine phosphatases. Available as a 1 ml vial or a set of five 1 ml vials. Once reconstituted, each vial contains the following components: A cocktail of four phosphatase inhibitors for broad-spectrum inhibition of acid, alkaline, serine/threonine and protein tyrosine phosphatases. Available as a 1 ml vial or a set of five 1 ml vials. Each vial, when reconstituted with 1 ml H<sub>2</sub>O contains 250 mM Sodium Fluoride, 50 mM beta-Glycerophosphate (Cat. No. 35675), 50 mM Sodium Pyrophosphate Decahydrate, 50 mM Sodium Orthovanadate.