

Tris Base, Molecular Biology Grade - CAS 77-86-1 - Calbiochem TRIS base is useful in the pH range of 7.0-9.0. Has a pKa of 8.1 at 25C.

Art. ID SAF-648310-2.5KG
Unit 1 x 2,5 kg
Deliverydetails No Dangerous Good /not restricted

Description

Useful in the pH range of 7.0-9.0. Has a pKa of 8.1 at 25C., Widely used as a biological buffer (absorbance: ≤ 0.03 at 260 nm, ≤ 0.02 at 280 nm). Useful in the pH range of 7.0-9.0. Has a pKa of 8.1 at 25C.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Trometamol	[77-86-1]				