

**PhosphoDetect Anti-Src Family (pTyr<sup>?</sup>) Rabbit pAb liquid, Calbiochem(R)**

Art. ID                      SAF-PK1109-50UL  
Unit                            1 x 50 µL

Description

Recognizes the ~60 kDa Src protein phosphorylated at Tyr416 in COLO 201 cells and PDGF-treated NIH/3T3 cells. May also recognize other Src family kinases phosphorylated at equivalent sites., This PhosphoDetect Anti-Src Family (pTyr<sup>?</sup>) Rabbit pAb is validated for use in Immunoblotting for the detection of Src Family (pTyr<sup>?</sup>), Protein A and immunoaffinity purified rabbit polyclonal antibody. Recognizes the ~60 kDa Src protein phosphorylated at Tyr416.