

PhosphoDetect Anti-JAK2 (pTyr1007/1008) Rabbit pAb liquid, Calbiochem(R)

Art. ID SAF-PS1014-100UL
Unit 1 x 100 µL

Description

This PhosphoDetect Anti-JAK2 (pTyr1007/1008) Rabbit pAb is validated for use in Immunoblotting for the detection of JAK2 (pTyr1007/1008)., Recognizes the ~130 kDa JAK2 protein phosphorylated at Tyr1007/1008 in LIF-treated 3T3-L1 cells., Rabbit polyclonal antibody adsorbed against non-phosphopeptide corresponding to the immunogen phosphorylation site, followed by immunoaffinity chromatography. Recognizes the ~130 kDa JAK2 protein.