

Custom Solution - 2 components 1000 mg/L each of in HNO3 5%

Art. ID LM24-CUS-71223-100ML
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		