

Custom solution - 13 components in 2% HNO₃

Art. ID LM24-CUS-71896
Unit 2 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	2	µg/mL		
	Zinc (Zn)	[7440-66-6]	20	µg/mL		
	Aluminium (Al)	[7429-90-5]	20	µg/mL		
	Iron (Fe)	[7439-89-6]	20	µg/mL		
	Vanadium (V)	[7440-62-2]	5	µg/mL		
	Chromium (Cr)	[7440-47-3]	5	µg/mL		
	Manganese (Mn)	[7439-96-5]	5	µg/mL		
	Nickel (Ni)	[7440-02-0]	5	µg/mL		
	Copper (Cu)	[7440-50-8]	5	µg/mL		
	Arsenic (As)	[7440-38-2]	1	µg/mL		
	Selenium (Se)	[7782-49-2]	1	µg/mL		
	Cadmium (Cd)	[7440-43-9]	1	µg/mL		
	Uranium (U)	[7440-61-1]	1	µg/mL		