

**SIB-Testblock - Vickers Hardness, DIN EN ISO 6507-3 and ASTM E 384, 65mmx10mm-HV,  
Procedure: HV1, Hardness: 300, Grid Lines**

|                 |  |
|-----------------|--|
| Art. ID         | SIB-VI-65x10-HV-6507-3/E384-HV1-300-GL |
| Unit            | each                                   |
| Deliverydetails | No Dangerous Good /not restricted      |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory