

**SIB-Testblock - Rockwell Hardness, DIN EN ISO 6508-3, 60mmx60mmx16mm-HR,  
Procedure: HRC, Hardness: 40, Grid Lines**

|                 |                                     |
|-----------------|-------------------------------------|
| Art. ID         | SIB-RW-60x60x16-HR-6508-3-HRC-40-GL |
| Unit            | each                                |
| Deliverydetails | No Dangerous Good /not restricted   |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory