

**SIB-Testblock - Rockwell Hardness, DIN EN ISO 6508-3 and ASTM E 18,  
60mmx60mmx16mm-HR, Procedure: HR45N, Hardness: 69, Grid Lines**

|                 |   |
|-----------------|---|
| Art. ID         | SIB-RW-60x60x16-HR-6508-3/E18-HR45N-69-GL |
| Unit            | each                                      |
| Deliverydetails | No Dangerous Good /not restricted         |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory