

**SIB-Testblock - Rockwell Hardness, DIN EN ISO 6508-3 and ASTM E 18,  
60mmx60mmx16mm-HR, Procedure: HR15N, Hardness: 87**

|                 |  |
|-----------------|--|
| Art. ID         | SIB-RW-60x60x16-HR-6508-3/E18-HR15N-87 |
| Unit            | each                                   |
| Deliverydetails | No Dangerous Good /not restricted      |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory