

**SIB-Testblock - Brinell Hardness, DIN EN ISO 6506-3, 100mmx100mmx16mm-HB,  
Procedure: HBW10/3000, Hardness: 300**

|                 |  |
|-----------------|--|
| Art. ID         | SIB-BR-100x100x16-HB-6506-3-HBW10/3000-300 |
| Unit            | each                                       |
| Deliverydetails | No Dangerous Good /not restricted          |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory