

**SIB-Testblock - Brinell Hardness, DIN EN ISO 6506-3 and ASTM E 10,  
100mmx100mmx16mm-HB, Procedure: HBW2.5/187.5, Hardness: 600, Grid Lines**

Art. ID                      SIB-BR-100x100x16-HB-6506-3/E10-HBW2.5/187.5-600-GL  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory