

**SIB-Testblock - Brinell Hardness, DIN EN ISO 6506-3, 100mmx100mmx16mm-HB,
Procedure: HBW2.5/62.5, Hardness: 250, Grid Lines**

| | |
|-----------------|--|
| Art. ID | SIB-BR-100x100x16-HB-6506-3-HBW2.5/62.5-250-GL |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory