

**SIB-Testblock - Brinell Hardness, DIN EN ISO 6506-3 and ASTM E 10,
100mmx100mmx16mm-HB, Procedure: HBW1/30, Hardness: 450, Grid Lines**

| | |
|-----------------|--|
| Art. ID | SIB-BR-100x100x16-HB-6506-3/E10-HBW1/30-450-GL |
| Unit | each |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Producer: DIN EN ISO 17025 DAkkS-accredited Laboratory