

Custom Solution - 5 components 10 µg/mL each of in Dichloromethane

Art. ID LM24-CUS-71919
Unit 1 x 5 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1593 III E1

Description

Conservation 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------|------------|---------------------|-------|--------|--------|
| | Benzo(b)fluoranthene | [205-99-2] | 10 | µg/mL | | |
| | Benzo(k)fluoranthene | [207-08-9] | 10 | µg/mL | | |
| | Benzo(a)pyrene | [50-32-8] | 10 | µg/mL | | |
| | Indeno(1,2,3-c,d)pyrene | [193-39-5] | 10 | µg/mL | | |
| | Benzo(g,h,i)perylene | [191-24-2] | 10 | µg/mL | | |