

Custom Solution - Formaldehyde CAS:50-00-0, 37 wt. % in H₂O, contains 10-15% Methanol as stabilizer (to prevent polymerization) EC:200-001-8

Art. ID LM24-CUS-71926
Unit 2 x 5 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3286 8 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Formaldehyde	[50-00-0]	37	%		